

## DTI processing for structural cortical networks analysis

# First use MRIConvert on dicoms to create *.nii*, *bvals* and *bvecs*.

# pre-process with eddy-correct which produces *data.nii.gz*

`eddy_correct <dti_scan.nii.gz> data 0`

# use brain extraction tool (bet) to extract brain from image and add `< -m >` to automatically

# create binary brain mask

`bet data.nii.gz nodif_brain -m`

# run dtifit with Fdt

**Fdt**

- in GUI select dtifit

# this should output the following files

# *dti\_V1* - 1st eigenvector

# *dti\_V2* - 2nd eigenvector

# *dti\_V3* - 3rd eigenvector

# *dti\_L1* - 1st eigenvalue

# *dti\_L2* - 2nd eigenvalue

# *dti\_L3* - 3rd eigenvalue

# *dti\_MD* - mean diffusivity

# *dti\_FA* - fractional anisotropy

# *dti\_MO* - mode of the anisotropy (oblate  $\sim -1$ ; isotropic  $\sim 0$ ; prolate  $\sim 1$ )

# *dti\_S0* - raw T2 signal with no diffusion weighting

# check results, e.g. and change DTI display to RGB and modulate with FA map

`fslview dti_V1 dti_FA`

# create forbedpost directory and copy data, nodif\_brain\_mask, nodif\_brain, bvals, bvecs

`mkdir forbedpost`

`cp -av bvals bvecs nodif_brain* data.nii.gz forbedpost/`

# run bedpostX and use forbedpost dir as input directory in GUI. This will create

# forbedpost.bedpostX directory. This will take approx 20 hours (up to a day).

**Fdt**

- Select BEDPOSTX

- input dir <forbedpost>

# this should create the following outputs in forbedpost.bedpostX

# *<i>* indicates the i-th fibre. It ranges from 1 to the maximum number of fibres set in the  
# advanced options.):

# **merged\_th<i>samples** - 4D volume - Samples from the distribution on theta

# **merged\_ph<i>samples** - Samples from the distribution on phi

# theta and phi together represent the principal diffusion direction in spherical polar

# coordinates

# **merged\_f<i>samples** - 4D volume - Samples from the distribution on anisotropic volume  
# fraction (see technical report).

# **mean\_th<i>samples** - 3D Volume - Mean of distribution on theta

**# mean\_ph<i>samples** - 3D Volume - Mean of distribution on phi  
**# mean\_f<i>samples** - 3D Volume - Mean of distribution on f anisotropy  
# Note that in each voxel, fibres are ordered according to a decreasing mean f-value  
**# mean\_dsamples** - 3D Volume - Mean of distribution on diffusivity d  
**# mean\_S0samples** - 3D Volume - Mean of distribution on S0 intensity  
**# dyads<i>** - Mean of PDD distribution in vector form. Note that this file can be loaded into  
# fslview for easy viewing of diffusion directions  
**# dyads<i>\_dispersion** - 3D Volume - Uncertainty on the estimated fibre orientation.  
# Characterizes how wide the orientation distribution is around the respective dyad.  
**# nodif\_brain** - brain extracted version of nodif - copied from input directory  
**# nodif\_brain\_mask** - binary mask created from nodif\_brain - copied from input directory

# change dir into forbedpost.bedpostX/  
# copy perl scripts to home/abcd1234/bin/ and add dir\_path (only first time)  
# aal\_images.pl generate\_matrix.pl calc\_total\_matrix.pl calculate\_vector.pl  
# calculate\_matrix.pl

# use flirt for getting nodif\_brain into standard space.  
# this will also create a nodif2std.mat file

[Flirt](#)

# in GUI use following  
- ref is standard set to T1\_2mm\_mni brain.  
- in nodif\_brain.nii.gz  
- out nodif2std

# invert the nodif2std.mat file

[InvertXFM](#)

# in GUI use following  
- in nodif2std.mat  
- out std2nodif.mat

# Apply inverse mat-file on aal.nii.gz to create aal image in dti space

[ApplyXFM](#)

# in GUI use following  
- mat std2nodif.mat  
- in /home/psyc0613/aal.nii.gz  
- ref nodif\_brain.nii.gz  
- out aal2nodif.nii.gz

**Use nearest neighbour !**

# Perform visual check on *data.nii.gz* and *aal2nodif.nii.gz*

# run first perl script to create 90 masks (only cortical areas from aal space).  
# this will also create a directory aal\_im where all 90 masks will be stored in the following  
# format: aal2nodif\_xx.nii.gz, furthermore a file called target\_list.txt will be created, which is  
# where the 90 masks are listed (check aal\_im dir for all 90 masks).  
# run the script from the forbedpost.bedpostX folder.

[aal\\_images.pl aal2nodif.nii.gz](#)

```
# generate connections with perl script. This will create 90 directories, 1 for each mask/region
# with all other regions in the form of XX_all
# this step will take several days
# the following files from bedpostX will be used:
# merged_th<i>samples
# merged_ph<i>samples
# merged_f<i>samples
# nodif_brain_mask
generate_matrix.pl

# check if *_all dirs all have 90 seeds_to_aal2nodif_XX.nii.gz, fdt_paths.nii.gz,
# probtrackx.log and waytotal,

# now run the ROI vector perl script. This will create a ROI_XX_probmatrix.txt in every
# XX_all dir. Output: Calculates the sizes of seed ROIs (number of voxels for each XX_all
dir.)
calculate_vector.pl

# Now generate a connectivity matrix, this will create a 90x90 matrix
calc_new_total_matrix.pl

# Now generate a probability matrix, this will create a different 90x90 matrix
calc_total_prob_matrix.pl
```